

In-house X-ray Standing Wave study of LaN/B multilayer mirrors

C.P. Hendriks, I.A. Makhotkin*, S.N. Yakunin**, R. van de Kruijs, F. Bijkerk


Industrial Focus Group XUV Optics, MESA+ Institute for Nanotechnology, University of Twente, Enschede, the Netherlands

** RNC Kurchatov Institute, Moscow, Russia

*E-mail: i.makhotkin@utwente.nl

Introduction

- X-ray Fluorescence (XRF) used for non-destructive characterization of atomic concentration of thin film.
- The X-ray standing wave (XSW) combines Grazing Incidence X-ray Reflectivity (GIXR) and XRF to reconstruct the atomic profile of thin films and multilayers.
- Measurements done on PANalytical Empyrean XRD setup equipped with Amptek XR-100 Silicon Drift Detector.



- Goal:
- Demonstrate applicability of XSW on in-house XRD setup.

- Benefits:
- Non-destructive
 - Suitable for ambient
 - Realizable on standard XRD equipment


The reconstruction of the multilayer period

- Samples: 2 periodic La/B and LaN/B multilayer mirrors:
- One baseline sample with clear intermixing between the 2 layers;
 - Passivation of La with nitrogen reduces intermixing between the La and B layers;
 - XSW was used to study the width of La distribution.


GIXR Measurements

The optical constant profiles are obtained from the GIXR data analysis;


- The passivated sample shows a more centered distribution.



Fluorescence Measurements




The fluorescence intensity is measured for different angles and the intensity per peak is obtained.



The angular fluorescence intensity shows a modulation by the standing wave

Electric Field Calculations


- Field Intensity is calculated using dynamical theory and the obtained optical profiles from GIXR data;
- The field is summarized to create an effective field for one effective period.
- A clear effective standing wave is visible in the first Bragg peak.



Calculation of Profiles

The profiles are calculated by solving the linear system (1);

- The passivated sample shows a more centered La distribution and more La is present in the B layer in the baseline sample;




→ The La barrier in the passivated sample has less intermixing with the B barrier than for the baseline sample → Expected result is a good confirmation of the effectiveness of the XSW technique.

Project: model independent approach to the GIXR data analysis

- Problem of classical of GIXR analysis: initial guess (typically bi-layer model) required
- Outcome of the fit depends on the initial model
- Time consuming manual procedure

Solution:

- Present a multilayer period as a set of "independent bins" and vary only the optical constants of these bins.
- No pre-defined number of sub-layers in ML period
- No pre-defined shape of the interfaces
- The same initial model for all analyzed multilayers
- Even useful in more complex systems.



Result: Fast and reliable fit!

Model independent approach to XSW data analysis

$I(\theta)$ -Measured fluorescence intensity:

$$I(\theta) = \sum_{j=1}^m P(z_j) |E(\theta, z_j)|^2$$

E - 2D Electromagnetic field

$P(z)$ -unknown atomic profile

To find $P(z)$: the system of linear equations to be solved:

$$\hat{A} \cdot \mathbf{x} = \mathbf{b} \quad (1)$$

Where: $\mathbf{x} = \begin{Bmatrix} P_1 \\ P_2 \\ \dots \\ P_m \end{Bmatrix}$


A = Electric field intensity for different angles and depths
 b = experimentally measured fluorescence

Measurements have errors and the problem is over-defined → ill-posed problem
 Solution → regularization required!


$$(\hat{A} + \lambda \cdot D) \cdot \mathbf{x} = \mathbf{b}$$

$$D = \begin{bmatrix} 2 & -1 & 0 & 0 & 0 & -1 \\ -1 & 2 & -1 & 0 & 0 & 0 \\ \dots & \dots & \dots & \dots & \dots & \dots \\ -1 & 0 & 0 & 0 & -1 & 2 \end{bmatrix}$$

λ - regularization coefficient



Experimental scheme: classical scheme for GIXR measurements + energy dispersive X-ray fluorescence detector.





Effective EM field from the periodic part of the La/B MLM in the first Bragg reflection condition (center graph), calculated GIXR (left graph) and fluorescence yield from La and B (right graph), and the δ profile used for calculation of EM field (top graph).

Application example: analysis of multilayer contamination

- XSW analysis of small Iron traces
- Iron can enter the system during deposition
- This can lead to a very small periodical contamination which is invisible to any non-atomic sensitive technique
- Sensitivity higher than a promille

Iron atoms in this multilayer were found in the interfaces.



Conclusion

- Model independent GIXR analysis allows a fast optical constant profile determination for following XSW analysis;
- The XSW technique is atomic sensitive which allows it also to differentiate between layers and materials that do not show any optical contrast;
- It is capable of discerning small differences in intermixing between 2 layers.
- Allows to find very small contaminations in thin films and also locate them.

1. Yakunin, S.N., et al., *Model independent x-ray standing wave analysis of periodic multilayer structures*. Journal of Applied Physics, 2014. **115**(13)
2. Makhotkin, I.A., et al., *Short period La/B and LaN/B multilayer mirrors for ~6.8 nm wavelength*. Optics Express, 2013. **21**(24)
3. Zameshin, A., et al., *Reconstruction of interfaces of periodic multilayers from X-ray reflectivity using a free-form approach*. Journal of Applied Crystallography, 2016. **49**(4)
4. A. N. Tikhonov, and I. A. Arsenin, "Solutions of ill-posed problems", Winston, (1977)

Acknowledgements. We acknowledge the support of the Industrial Focus Group XUV Optics at the MESA+ Institute at the University of Twente, notably the industrial partners ASML, Carl Zeiss SMT AG, PANalytical, SolMates, TNO, and Demcon, as well as the Province of Overijssel and the Foundation FOM.